

FACS-array protocol

Striatal Dissociation Protocol for FACS

Dissection

Set up equipment for dissection and vibratome for cutting (make sure to pack vibratome with ice to keep brain cold)

1. Anesthetize animal and dissect out brain (use 3-4 EGFP mice and 1 wt mouse)
2. Immediately place brains in slice solution (bubbled with carbogen)
3. Place brains on petri dish(w/ filter paper) and cut off cerebellum and tip of forebrain with razor blade
4. Crazy glue brain to vibratome, cover in slice solution and section at 400 μ m (EGFP brains can be cut together)
5. Remove sections to Petri dish (keep in slice solution) and dissect out striatum under dissecting microscope
6. Place striatal sections in papain mix

Dissociation

Prior to dissection:

1. If using a new albumin ovomucoid inhibitor (AOI) vial then add 32 ml EBSS and equilibrate with carbogen
2. Add 5 ml EBSS to papain mix and split mix into separate tubes for wt and EGFP tissue
3. Incubate papain mix at 37 $^{\circ}$ C while equilibrating with carbogen (95% O₂/5% CO₂)
4. Add 500 μ l EBSS to DNase vial and mix gently (keep DNase on ice until ready to use)

After dissection:

1. Add DNase to papain mix (50 μ l per 1 ml papain mix) **DNase is added to prevent the cells from clumping together**
2. Place striatal sections in papain mix and displace air in vial with carbogen
3. Incubate at 37 $^{\circ}$ C for 40min
4. Periodically shake tissue sections by hand (ideally a rocker platform should be used)
5. Triturate with 3 glass pipettes of decreasing tip diameter
6. Centrifuge @ 1200rpm for 5min
7. Prepare EBSS/AOI solution: 2.7ml EBSS , 300 μ l reconstituted AOI, and 150 μ l DNase
8. Remove supernatant and resuspend cell pellets in EBSS/AOI solution
9. Prepare discontinuous density gradient: Add 2-3ml AOI to a 15 ml tube and carefully layer cell suspension atop
10. Centrifuge @ 900rpm for 6min

11. Discard supernatant and resuspend pellet in 2ml FACS Buffer
 12. Pass cell suspension through 70µm mesh and collect in a 50 ml tube
 13. Centrifuge at 1200rpm for 5 min
 14. Remove supernatant and resuspend in 1 ml FACS Buffer
 15. Transfer to sterile FACS tube and add 2µl PI (propidium iodide) (1mg/ml) to stain dead cells red
 16. Transfer samples on ice to FACS core for sorting
- *keep cells on ice in between centrifugation and use a 4°C centrifuge if possible**

note: addition of phenol red can interfere with the lasers' ability to properly detect fluorescence

FACS Buffer recipe (50ml)

- 50mg BSA
- 0.5ml 1M Hepes
- 0.5ml 100x Pen-Strep
- 1.25mg DNase (2,020 units/mg)
- Bring up to 50ml with L15-CO2 w/ no Phenol Red

Slice Solution

| <u>Agent</u> | <u>mM</u> | <u>MW (mols)</u> | <u>2L</u> | <u>4L</u> |
|----------------------|-----------|------------------|-----------|-----------|
| NaHCO3 | 26 | 84.01 | 4.3686g | 8.7372g |
| NaH2PO4 . H2O | 1.25 | 137.99 | .3450g | .6900g |
| Glucose | 10 | 180.23 | 6.040g | 7.2080g |
| 5M NaCL | 130 | 58.44 | 52mL | 104mL |
| 2M KCl | 3 | 74.56 | 3mL | 6mL |
| 1 M MgCl2 | 5 | 95.21 | 10mL | 20mL |
| 1M CaCl2 | 1 | 147.02 | 2mL | 4mL |

ph=7.2-7.4
Osmo=290-310

Papain Dissociation Kit Vials

- Vial 1: EBSS
- Vial 2: Papainwith L-cysteine and EDTA
- Vial 3:DNase
- Vial 4: AOI

FACS/RNA extraction

1. FACS cells immediately- prolonging FACS has resulted in decreased EGFP+ cell recovery

2. We use the UCLA Flow Cytometry Core: There are multiple technicians sorting on multiple machines but it is best to use the same technician/machine for consistency
3. Use wt (non-EGFP) cells to set parameters
4. FACS set up should be set to detect PE (PI signal) and FITC (EGFP signal)
5. Keep cells at 4°C during FACS
6. Collect cells in FACS buffer at 4°C

We collect two populations of cells: Low EGFP expression- FITC signal 10^3 - $10^{3.5}$, High EGFP expression- FITC signal $>10^{3.5}$ (Total High EGFP recovery is typically 5,000-12,000 cells)

7. Spin down cells immediately following FACS and proceed with Arcturus Picopure RNA extraction protocol (w/ Qiagen DNase treatment)
8. Check RNA integrity on Bioanalyzer picochip (Geschwind Lab)

We use high EGFP expressing cells for further study- this population shows consistent RT-PCR results for our known genes.

Two Round RNA Amplification

Amplify and label with Agilent Low RNA Input Fluorescent Linear Amplification Kit (starting RNA concentration is approximately **3ng**- a relative concentration was obtained from the picochip. 2 rounds of amplification is necessary.)

* *For RT-PCR one round of amplification is sufficient.

DAY I

cDNA Synthesis from Total RNA

| Annealing Rxn | uL |
|--------------------|-------|
| RNA soln (3ng) | 10.3 |
| T7 Promoter Primer | 1.2 |
| <hr/> | <hr/> |
| Total Vol | 11.5 |

Incubate 65C 10 min, ice 5 min.

| cDNA Rxn | uL |
|------------------------|-------|
| Annealing Rxn | 11.5 |
| 5X First Strand Buffer | 4 |
| 0.1 M DTT | 2 |
| dNTP Mix (10mM) | 1 |
| MMLV RT | 1 |
| RnNase Out | 0.5 |
| <hr/> | <hr/> |
| Total Vol | 20 |

Incubate 40C
2hrs.

Inactivate 65C 15 min., ice 5 min.

IVT Amplification - Round 1



IVT Amplification Rxn

| | |
|---------------------|------|
| cDNA rxn | 20 |
| Nuclease free water | 12.1 |
| 4X Transcription | |
| Buffer | 20 |
| 0.1 M DTT | 6 |
| NTP Mix | 8 |
| CTP | 5.6 |
| 50% PEG | 6.4 |
| RnNase Out | 0.5 |
| Inorganic Phosphate | 0.6 |
| T7 RNA Polymerase | 0.8 |
| <hr/> | |
| Total Vol | 80 |

Incubate 40C for 16 hours (**Protocol suggested 4 hours to overnight. I did 16 hours the first time and used the same time for all my samples.**)

Purification by Qiagen RNAeasy Mini Kit, Elute in 30 uL Water
Vacuum centrifuge down to 10.5uL

(Recently I've been using Qiagen Rneasy mini elute kit to avoid the vacuum centrifuge.)

DAY II

cDNA Synthesis from cRNA

First Strand Synthesis

Annealing Rxn

| | |
|-----------------|------|
| cRNA Soln | 10.5 |
| Random Hexamers | 1 |
| <hr/> | |
| Total Vol | 11.5 |

Incubate 65C 10 min, ice 5 min.

cDNA Rxn - First Strand

| | |
|------------------------|------|
| Annealing Rxn | 11.5 |
| 5X First Strand Buffer | 4 |
| 0.1 M DTT | 2 |
| dNTP Mix (10mM) | 1 |
| MMLV RT | 1 |
| RnNase Out | 0.5 |
| <hr/> | |
| Total Vol | 20 |

Incubate 40C for 2 hrs.

Second Strand Synthesis

Annealing Rxn

| | |
|--------------------|------|
| cDNA Soln | 20 |
| T7 Promoter Primer | 1.2 |
| <hr/> | |
| Total Vol | 21.2 |


Incubate 65C 10 min, ice 5 min.

cDNA Rxn - Second Strand

| | |
|---------------|------|
| Annealing Rxn | 21.2 |
|---------------|------|

| | |
|------------------------|-----|
| Nuclease free water | 1.8 |
| 5X First Strand Buffer | 8 |
| 0.1 M DTT | 4 |
| dNTP Mix (10mM) | 2 |
| MMLV RT | 2 |
| RnNase Out | 1 |
| <hr/> | |
| Total Vol | 40 |

Incubate 40C
2hrs.
Inactivate 65C 15 min., ice 5 min.



****Some of this sample can be used for RT-PCR validation**

DAY III

IVT Amplification - Round 2

| IVT Amplification Rxn | Cy3 | Cy5 |
|------------------------------|------------|------------|
| cDNA rxn | 20 | 20 |
| CTP-Cy3 | 2.4 | 2.4 |
| Nuclease free water | 15.3 | 15.3 |
| 4X Transcription Buffer | 20 | 20 |
| 0.1 M DTT | 6 | 6 |
| NTP Mix | 8 | 8 |
| 50% PEG | 6.4 | 6.4 |
| RnNase Out | 0.5 | 0.5 |
| Inorganic Phosphate | 0.6 | 0.6 |
| T7 RNA Polymerase | 0.8 | 0.8 |
| <hr/> | | |
| Total Vol | 80 | 80 |

Incubate 40C for 5 hours

(Protocol suggested 4 hours to overnight. I did 5 hours the first time and used the same time for all my samples.)

Purification by Qiagen RNAeasy Mini Kit, Elute in 60 uL (2 X 30 uL) Water

- Prior to hybridization check concentration on Nanodrop and Bioanalyzer nanochip (Geschwind lab)
- Hybridize to mouse and mouse developmental Agilent arrays (Microarray core)

FACS-array reagents/kits/etc

| Product | cat# | company | price (\$) |
|--|-------------|---------------------|-------------------|
| Cell Dissociation/FACS Dissociation | | | |
| Papain Dissociation System (1bx) Slice solution (see protocol for recipe) | LK003150 | Worthington Biochem | 160.00 |

glass pipettes

FACS Buffer reagents

| | | | |
|---------------------------------|-----------|------------|--|
| Leibovitz-L15 media (no phenol) | 21083-027 | Invitrogen | |
| BSA | | | |
| Pen-Strep 100x | P0781 | Sigma | |
| Dnase I | D 4527 | Sigma | |

Hepes

Other reagents

| | | | |
|---------------------------------|-----------|-----|------------|
| Propidium Iodide | | | |
| 70µm mesh cell strainer | 21008-952 | VWR | |
| FACS (UCLA Flow cytometry core) | | | \$70/hour? |

RNA

RNA extraction

| | | | |
|-------------------------|---------|----------|--|
| Picopure RNA Extraction | KIT0204 | Arcturus | |
|-------------------------|---------|----------|--|

RNA quality

| | | | |
|--------------------------------|-----------|---------|--------|
| Picochip (Agilent Bioanalyzer) | 5065-4473 | Agilent | 598.00 |
| Nanochip (Agilent Bioanalyzer) | 5065-4476 | Agilent | 487.00 |

RNA Amplification

| | | | |
|------------------------------------|-----------|------------|--------|
| Low Input Linear Amp kit (20 rxns) | 5184-3523 | Agilent | 798.00 |
| Cy-3 CTP (100nmol) | NEL580 | Perkin Elm | 214.20 |
| Cy-5 CTP (100nmol) | NEL581 | Perkin Elm | 214.20 |
| Rneasy Mini kit | 74104 | Qiagen | 203.00 |

(alternative : order Rneasy mini elute kit so RNA doesn't have to be vacuum centrifuged after 1st round)

Microarray Hybridization

| | | | |
|---|-------------|---------|---------|
| Agilent Whole mouse oligo Arrays (5 slides) | G4122A | Agilent | 2875.00 |
| Hybridization kit Plus (for 10 slides) | 5184-3568 | Agilent | 226.00 |
| Stabilization and Drying Solution, 500 mL | 5185-5979 | Agilent | 95.00 |
| Surehyb backings (5 slides) | G2534-60003 | Agilent | 100.00 |
| Other | | | |
| Scanning (per array) | Array Core | Agilent | 25.00 |